


GEIGER UNIVERSAL MOBILE STAND: PART NO. 150-MS

CAUTION: FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN OR HEALTH CARE PROFESSIONAL.



Inspection

Inspect the Universal Mobile Stand upon receipt. If damage has occurred during shipment, do not use. Notify Geiger immediately.

Package should contain: (1) base w/5 casters, (2 locking), (1) pole with attached back plate.

Tools Needed:

1/8 inch Allen wrench

Assembly Instructions

1. Position post (A) into hub of base.
2. Insert pole into base, tighten both set screws until tight.
3. Note: The Universal Mobile Stand requires the 150-W Wall Mount Kit (sold separately) in order to attach the device to the mobile stand.

